Fiber Optic Cables for UV-VIS-NIR

- Single strand, multimode fiber optic cables are available with SMA905 connectors, ¼ inch dia. steel mono-coil armor with black PVC jacket. Sizes include 400, 600, or 1000μm core diameters
- > Standard lengths are 1 and 2 meters
- > Custom lengths are available with a minimum price of a 1 meter fiber (extra delivery time also)
- > Solarization resistant fibers available for deep UV applications below 300nm in all diameters
- ➤ **High OH fibers** (190-850 nm) are needed for UV and UV-VIS applications below 400nm
- ➤ Low OH fibers (400-2200 nm) are needed for VIS-NIR applications. They offer lowest possible internal light attenuation. Intended for applications with long fiber lengths.

> Specialty fibers for various applications

- Y-cables to couple 2 spectrometer systems
- Y-bundled cables to couple 2 light sources
- 1 to 3 splitter cable
- Custom fibers to meet your application



Armored Fiber Optic Cable Splitter 1-TO-3

400 µm Fibers

| Item | Description | Price |
|----------------|---|-------|
| F400-UVVis-SR | Armored, 2 meters, Solarization Resistant fiber optic cable | \$175 |
| F400-VISNIR | Armored, 2 meters, fiber optic cable | \$155 |
| F400-per meter | Price per additional meter, armored with black PVC jacket | \$70 |
| F400YBNIRUVSR | Y-cable bundle, 0.5 meter, 4 SR & 3 NIR fibers connects SL1+SL3 lamps | \$380 |

600 µm Fibers

| Item | Description | Price |
|----------------|---|-------|
| F600-UVVis-SR | Armored, 2 meters, Solarization Resistant fiber optic cable | \$190 |
| F600-VisNIR | Armored, 2 meters, fiber optic cable | \$160 |
| F600-per meter | Price per additional meter, armored with black PVC jacket | \$110 |
| F600-Y-VISNIR | Armored Y cable, 1 meter, 1 SMA to 2 SMA, all for NIR | \$280 |
| F600-Y-UVSRNIR | Armored Y cable, 1 meter, 1 SMA to 2 SMA, SR+NIR | \$380 |
| F600-1TO3-SR | Armored splitter, 1 meter, 1 SMA to 3 SMA, all Solarization Resistant | \$365 |

1000 µm Fibers

| Item | Description | Price |
|------------------|--|-------|
| F1000-UVVis-SR-1 | Armored, 1 meter, 1000um, Solarization Resistant fiber optic cable | \$210 |
| F1000-VISNIR-1 | Armored, 1 meter, 1000um fiber optic cable | \$180 |
| F1000-per meter | Price per additional meter, armored with black PVC jacket | \$150 |

| Specifications | | | |
|-----------------------|---------------------------|-------------------|-----------------------------|
| Weight | Varies on length/diameter | Fiber material | Silica core - silica clad |
| Length | 1-2 meter, or custom | Diameter | 400 μm or 600 um or 1000 μm |
| Spectral range | 190 – 2200nm | Connectors | SMA 905 |
| Numerical aperture | 0.22 NA | Temperature range | 0 – 50 Deg C |