

... ONE APP TO TREAT THEM ALL ...







Spectragryph 1.2 software for optical spectroscopy

Tired & Exhausted of:

- Booking spectrometer time just to process your spectral data?
- Stupefying conversion of original spectral files all the time?
- Having to use Origin or even Excel to bring your spectra together?

Put an end to your suffering with the new Spectragryph software!

- 63 file formats (*Thermo, Jcamp-DX, Omnic, Bruker, Jasco, Shimadzu, Andor, PerkinElmer, Agilent, Horiba, Avantes, OceanOptics, StellarNet, B&W Tek, ...*)
- Ready for UV-VIS, NIR, FTIR, Raman, LIBS, XRF, fluorescence data
- Easy to use. Rock solid. Lean. Unbeatable. Powerful. Instant plots. Batch export.
- Process multiple files with automated sequences. Spectral libraries. Acquisition.



Spectroscopy Ninja – providing optical spectroscopy software since 2002 (*free* for academic use)

<http://spectroscopy.ninja> © 2016-2018: Dr. Friedrich Menges Software-Entwicklung